

# The EleGARD™ Patient Positioning System

The EleGARD™ Patient Positioning System (EleGARD) is a cardiopulmonary board which may elevate a patient's head and thorax: including during airway management; during manual CPR, manual CPR adjuncts, CPR with the LUCAS® Chest Compression System or ARM XR Automated Chest Compressor; and patient transport. EleGARD is indicated for adults only when used with the LUCAS Chest Compression System or the ARM XR Automated Chest Compressor. When not used with an automated compression system, EleGARD can be used for adults and children, not including infants and neonates.



List price of \$8,495 includes:

- 1 EleGARD Device
- 1 LUCAS 3 Backplate
- 2 Rechargeable Lithion Ion Batteries
- 1 Battery Charger
- 1 Carrying Case
- Instructions for Use

Disposable covers are available for purchase in packs of 10 for list price \$99.50

Weight: ~20 pounds with LUCAS Backplate and battery installed.

Size: ~30" long, ~11" tall, ~19" wide (in lowered position for storage.)

**For additional information or to contact a representative: 763.259.3722 or [info@elevatedcpr.com](mailto:info@elevatedcpr.com)**



**ElevatedCPR.com**

US and Foreign Patents: [elevatedcpr.com/patents](http://elevatedcpr.com/patents)  
©2023 – AdvancedCPR Solutions®

Maintenance and Service information: Apart from cleaning, disinfection and ensuring the EleGARD System battery is fully charged, the product does not require routine preventive maintenance or periodic calibration by a qualified technician. Out of warranty repair service is available from an ACS qualified service partner upon request. Battery packs, battery chargers, carrying cases, and neck foam replacement parts are available for purchase, and replacement part pricing is available upon request.